

Faces Patterns & Vivid Colors Palettes



U.S. Flag
66001-000/1



Water
66002-000/1



Candy Corn
66003-000/1



Camouflage
66004-000/1



Psychedelic
66005-000/1



Leopard
66006-000/1



Rainbow
66007-000/1



Flames
66008-000/1



Marble
66009-000/1



Red
66010-000/1



Yellow
66011-000/1



Green
66012-000/1



Blue
66013-000/1



Orange
66014-000/1



Pink
66015-000/1



Purple
66016-000/1




Black
66017-000/1



White
66018-000/1

Faces




Wanna be cool?
Snap on Vivid Blue.



Hear The
Call of The Wild?
Switch to Leopard.



Got a Hanking
for Sweets?
Try Candy Corn.



Feeling Groovy?
Go Psychedelic!

Faces

Put a Face On Better Hearing!

Now better hearing can be hot or cool, outrageous or discreet, camouflaged or flamboyant. All with Faces - the unique new interchangeable coverings for J13 behind-the-ear (BTE) hearing instruments.

With their fast, snap-on design, Faces let you "customize" your BTE's with one of 36 colors and patterns. Ideal for pediatric fittings, they can help kids develop unique and positive attitudes about amplification. And they provide a wide range of hair and skin tones for appearance-conscious patients of any age.



Faces' easy snap-on design lets you quickly change your BTE's outer appearance!

Faces offer much more than an expression of your individuality.

They can enhance instrument durability and moisture resistance by providing a protective covering over three sides.

Faces are available in two configurations. One allows access to all controls, the other covers everything but the microphone opening.

Available individually and in nine-pair display kits, Faces is the unique, convenient way to heighten or reduce BTE visibility.



All of Faces' colors and patterns are available in two opening configurations.



Faces' 9 pair display kit

Faces Human Color Palette

Human Color #1 66019-000/1	Human Color #2 66020-000/1	Human Color #3 66021-000/1
Human Color #4 66022-000/1	Human Color #5 66023-000/1	Human Color #6 66024-000/1
Human Color #7 66025-000/1	Human Color #8 66026-000/1	Human Color #9 66027-000/1
Human Color #10 66028-000/1	Human Color #11 66029-000/1	Human Color #12 66030-000/1
Human Color #13 66031-000/1	Human Color #14 66032-000/1	Human Color #15 66033-000/1
Human Color #16 66034-000/1	Human Color #17 66035-000/1	Human Color #18 66036-000/1