

We're pleased to offer AristaLite digital hearing technology as one of our many effective hearing solutions.

Don't let hearing loss keep you from hearing life's important sounds. Ask us today for more information about AristaLite.

To learn more about your hearing and available options visit [www.starkey.com](http://www.starkey.com).

# ARISTA *lite*™

## A Perfect Balance of Natural and Amplified Sound.

**Your Authorized Professional:**




## Why wear hearing instruments that amplify sounds that you already hear normally?

The new AristaLite's ideal blend of digital signal processing and open-ear acoustics delivers high frequency amplification without altering the lower frequency sounds you have no problem hearing.

The result is a more natural balance and clarity of sounds with greater listening comfort. Plus, common frustrations like "my voice sounds hollow" or feeling "plugged up" are greatly reduced.



Easily concealed behind the ear, AristaLite's ultra-sleek design also incorporates Starkey's exclusive Precision Audio Imaging. With its seven-band equalizer, your hearing professional can fine-tune AristaLite's high fidelity circuit to your exact listening needs. This reduces the risk of over-amplification and distortion – a problem particular to many open-ear fittings.



AristaLite's natural sound quality is further enhanced by Precision Audio Imaging's expansion and anti-feedback capabilities. These flexible algorithms help manage the sensitivity of digital signal processing in specific acoustic environments.

For added flexibility, AristaLite lets you switch between three standard



memories. With a simple button push, you can adjust the programmable settings to provide ideal clarity in a variety of sound environments.

Ask your hearing professional today about how AristaLite's unique blend of natural and amplified sound might provide you with superior high frequency hearing.