



Nine companies. One wireless solution.

In a unique cooperation, nine industry leaders – Audibel, Bernafon, Micro-Tech, Nu-Ear, Omni, Oticon, Phonic Ear, Qualitone, and Starkey – are creating a new frontier in personal wireless communication with Lexis™.


Freedom in every situation

Lexis operates with three levels of directionality and can be handheld, worn around the neck, placed on a tabletop, or connected to entertainment equipment.

Improved speech understanding

Combined with BTE hearing instruments, Lexis helps users in listening situations where obstacles like noise, distance, and reverberation make speech understanding difficult.

The Lexis microphone combines a directional four microphone array with advanced digital signal processing to provide outstanding directionality of up to 8.5 dB AI-DI. This allows the user to filter speech from unwanted noise, resulting in optimal speech understanding.



three listening modes

Simple, flexible fitting

A hearing care professional can easily connect and fit the ear-level Lexis FM receiver to virtually any hearing aid via an audio shoe. With Lexis it is also possible to adjust receiver sensitivity on-site on an ear-level receiver.

Intelligent design

The transmitter/microphone can be recharged like a mobile phone, but alkaline batteries may also be used. Attractive and ergonomically sound, the transmitter/microphone was designed with the end-user in mind.

LEXIS™

Nine companies have come together to bring Lexis to hearing care professionals and the patients they serve.

Lexis is the first in a line of inoperable, convenient products for better speech intelligibility and greater participation in life.

For more information, contact your preferred supplier:

