

# Understand the most important conversations with PDI Hearing Instrument Technology

## CONVENTIONAL INSTRUMENTS

Conventional aids amplify sounds from all directions...



...making it more difficult to understand a conversation with someone facing you.

## PDI ENHANCED DIGITAL INSTRUMENTS

When activated, PDI narrows the sound reception area...



...reduces interference of sounds and conversations from the sides and rear.




## IMPROVED HEARING WITH PDI


PDI's digital signal processor continuously monitors your sound environment. When noise levels reach a higher volume, the built-in computer narrows the sound pattern to focus on the area in front of you. Understanding in noisy places can be significantly improved, especially compared to hearing aids without this directional capability.



[www.starkey.com](http://www.starkey.com)



**PDI ITC**  
from Starkey Labs



**PDI ITE**  
from Starkey Labs



**PDI BTE**  
from Starkey Labs

# DEALER IMPRINT AREA