

Understand the most important conversations with PDI Hearing Technology



CONVENTIONAL INSTRUMENTS

Conventional aids amplify sounds from all directions...



...making it more difficult to understand a conversation with someone facing you.

PDI ENHANCED DIGITAL INSTRUMENTS


When activated, PDI narrows the sound reception area...




...reduces interference of sounds and conversations from the sides and rear.

IMPROVED HEARING WITH PDI FROM STARKEY

PDI's digital signal processor continuously monitors your sound environment. When noise levels reach a higher volume, the built-in computer narrows the sound pattern to focus on the area in front of you. With fewer sounds competing for your attention, understanding in noisy places can be significantly improved, especially compared to hearing aids without PDI.



PDI ITC



PDI ITE



PDI BTE



www.starkey.com

DEALER IMPRINT AREA