

Hear Better in Noise with Precision Directional Imaging™ (PDI)



VS



CONVENTIONAL DIGITAL HEARING INSTRUMENTS

Conventional aids amplify sounds from all directions making it more difficult to understand a conversation with someone facing you.

STARKEY'S PDI ENHANCED DIGITAL INSTRUMENTS

PDI narrows the sound focus area and reduces interference of competing sounds and conversations from the sides and rear.

IMPROVED HEARING WITH PDI FROM STARKEY

PDI's digital signal processor continuously monitors your sound environment. When noise levels reach a higher volume, the built-in computer narrows the sound pattern to focus on the area in front of you. With fewer sounds competing for your attention, understanding in noisy places can be significantly improved, especially compared to hearing aids without PDI.



PDI ITC



PDI ITE



PDI BTE



www.starkey.com

DEALER IMPRINT AREA