

Instead of hearing ALL the conversations...




...understand the one that matters most.

Restaurants and parties are full of sounds and conversations you don't necessarily want to hear. That's why we offer **PDI - Precision Directional Imaging**. An elaborate digital

processing system that narrows your hearing instruments' sound reception patterns to focus on the area directly in front of you. And the conversation that's most important.



CONVENTIONAL AMPLIFICATION



PRECISION DIRECTIONAL IMAGING



www.starkey.com

For further hearing information, ask your hearing professional and visit our website