

N - HI LO LIGHT PHOTO ON ON ON ON

WELCOME TO
GENESIS 4

AXENT 4

ARISTA 3

SEQUEL 2

ENDEAVOUR 1

MOVE YOUR PATIENTS UP TO THE NEXT LEVEL OF DIGITAL PERFORMANCE.

Starkey is proud to introduce **GENESIS 4**. Our family of new and enhanced digital signal processing solutions featuring Precise Audio Imaging for every level of patient listening and lifestyle need.

At the top, there's **AXENT II**, our premier four channel, eight band digital technology. Axent II features dynamic directionality and expanded feedback cancellation, with multiple setting noise management and improved expansion capabilities.

Also new is **ARISTA**, our three channel system with seven band, compression-independent filter shaping. Arista also offers directionality, noise and feedback management, as well as expansion.

SEQUEL, our two-channel, and **ENDEAVOUR**, our entry-level single channel technologies both now feature low-level

expansion as a part of their popular performance and efficiency features.

GENESIS 4 products also offer varying levels of convenience features like in-situ audiometry and verification, indicator tones for low battery, memory, and volume settings, as well as the industry's highest levels of power matrices. And all four utilize our proven compression architecture and unique high fidelity, yielding maximum audibility and comfort.

These products' capabilities are complemented by **PFS version 4.0**, Starkey's versatile Professional Fitting System, with its many productivity and performance enhancements.

Let **GENESIS 4** bring you and your patients to whole new levels of digital hearing excellence.

INTRODUCING GENESIS 4

STARKEY LABORATORIES INC. • www.starkey.com • 800.328.8602 • 952.941.6401

© ALL RIGHTS RESERVED • TJAD7256-00-EE-ST



VISIT US AT AAA 2003
TO SEE OUR NEWEST
TECHNOLOGY & SERVICES