



Information below is for all technology levels of Wi Series products: i110, i90 and i70.

Battery Life Expectancy Chart			
Standard Products	Life Expectancy (in days)*		Batteries Required (in 12 months)
	Without Streaming	With Streaming	
<b>Size 312</b>			
Wi Series RIC 40	9	6	41 - 61
Wi Series RIC 50	9	6	41 - 61
Wi Series RIC AP 60	9	6	41 - 61
Wi Series RIC AP 71	7	4	52 - 91

\*Assumes average current drain is 4.5 mA when streaming for 40, 50, 60 gain and 9.0 mA for 71 gain.  
All estimates are based on a 16-hour day and assumes 2 of those hours are spent streaming