



SOUNDGEAR

**Arm yourself
with the best**





Trust your ears to SoundGear

Whether you're in the field, at the range or on the factory floor, hearing is key to your success and safety. Yet the competing need to protect your hearing from the damaging sound of gun shots, power tools and industrial machines creates a unique noise management challenge.

SoundGear solves that challenge.



High-definition sound enhancement & automatic noise suppression

SoundGear represents the latest advancements in hearing protection. This unique family of products protects your hearing from loud noises like gun shots, engines and power tools and our electronic products enhance environmental awareness.

Made by Starkey — America's premier provider of hearing solutions — SoundGear devices feature advanced digital technology capable of enhancing sound eight times over traditional hearing protection products. That same technology simultaneously detects loud impulse noises as they happen, instantly and automatically suppressing them to protect your hearing.

SoundGear Features:

Dynamic digital noise protection for your ears

High-definition sound enhancement

Hear and locate flushes, whistling wings, gobbles and movement while also suppressing gun blasts

Noise Reduction Rating (NRR) of up to 27 dB*

Seamless sound-activated compression automatically suppresses noises

All-day, easy-to-wear comfort

*depending on model



PHANTOM

Introducing the world's first custom, Bluetooth® compatible and rechargeable hearing aid and protection device.



Phantom delivers benefits no other single product can:

- Pre-Programmed to the most common hearing loss range (Normal to Moderate)
- Automatically protect your ears from loud sounds and gun blasts
- Hear conversations, game and environmental sounds with ease
- Custom molded to your ear for the most precise and personalized fit
- Stream calls, music or other audio wirelessly from your phone
- One, overnight charge ensures you are powered all day long



Our complete kit contains:

- 1 pair of SoundGear Phantom devices
- 1 power adapter
- 1 charging case
- Zipper carrying case
- SoundGear removable lanyard system
- 1 cleaning brush
- 2 packs HearClear wax guards
- 2 SoundGear stickers

 **One-year, Worry-Free repair warranty**

 **30-day, risk-free trial***

TECHNICAL SPECIFICATIONS

NRR	Amplification	Suppression at
22 dB	30 dB	90 dB

- Volume Control
- 4 Digital Listening Modes (Normal, Max Compression, High Frequency Boost, Mute)
- Maximum Wind Noise Management applied to all programs
- For IPIL data see Phantom product page on SoundGear.com

*30-day, risk-free trial starts from receipt date

Electronic Custom Fit



SoundGear custom products are for those who refuse to compromise when it comes to performance and personalization. Each pair is custom molded to your ear.

Our complete kit contains:

- 1 pair of SoundGear custom electronic hearing protection devices
- 2 packs of batteries (Size 13)
- Zipper carrying case
- Black removable lanyard system
- 1 cleaning brush
- 2 packs HearClear wax guards
- 2 SoundGear stickers



One-year, Worry-Free repair warranty



30-day, risk-free trial*

TECHNICAL SPECIFICATIONS		
HUNTER & SHOOTER MODEL		
NRR	Amplification	Suppression at
24 dB	30 dB	95 dB
SILVER	PLATINUM	
• Volume Control	• Volume Control	
• 1 Digital Listening Mode (Normal)	• 4 Digital Listening Modes (Normal, High Frequency Boost + Max Wind Management, Telecoil, Mute)	

*30-day, risk-free trial starts from receipt date

Electronic Instant Fit In-the-Canal



SoundGear In-the-Canal is the smallest and lightest dynamic digital hearing protection product on the market. It's ideal for the hunter or shooter looking for an edge — or industrial workers looking for all-day comfort. Ready to wear right out of the box, it rests discreetly inside your ear to deliver natural wind reduction and superior sound quality.

Our complete kit contains:

- 1 pair of SoundGear electronic hearing protection devices
- 2 pairs of orange silicone sleeves – (1) small (1) large
- 2 packs of batteries (Size 10)
- Zipper carrying case
- 1 cleaning brush
- 1 SoundGear sticker



One-year, Worry-Free repair warranty



30-day, risk-free trial*

TECHNICAL SPECIFICATIONS		
HUNTER & SHOOTER MODEL		
NRR	Amplification	Suppression at
25 dB	15 dB	93 dB
INDUSTRIAL MODEL		
NRR	Amplification	Suppression at
25 dB	8 dB	80 dB

Non-Electronic Custom Fit

As the recognized experts in custom hearing instruments, Starkey and SoundGear integrate the same breakthrough manufacturing and technology innovation into our earplug solutions that are customized for every patient. **The result is unrivaled fit and comfort.**



Solid Ear Plugs

Custom-fit solid earplugs are made with super soft, high strength silicone and provide outstanding noise protection in a variety of recreational and occupational environments where high levels of noise are present.



Filtered Ear Plugs

Commonly referred to as musician or industrial earplugs, these custom-fit earplugs attenuate across all frequencies rather than just low- and mid-range frequencies, making music and environmental sounds much clearer and more natural. These are available with multiple filter options including a 10 dB, 17 dB, 26 dB, Impulse (33.1 dB IPIL) or solid filter.



High-Frequency Noise Filters (Hocks Noise Braker)

These custom earplugs contain a patented acoustical chamber and filter to attenuate sound to appropriate levels in the work place or on the range.





SOUNDGEAR

Global Headquarters
6700 Washington Ave. S.
Eden Prairie, MN 55344
1.800.769.8605

www.SoundGear.com

©2021 Starkey Laboratories, Inc. All Rights Reserved. 81372-007 7/21
BROC2663-05-EE-XX

Starkey, SoundGear, SoundGear logo and Phantom are registered trademarks of Starkey Laboratories, Inc.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks is under license.